

ANIMAL DATA

NEUROCAP[®]

Effective Nerve-End Barrier



POLYGANICS

TRANSFORMING PATIENT
RECOVERY

www.polyganics.com

NEUROCAP[®]

The Only Dedicated Nerve-end Capping Device

NEUROCAP[®]

Management of symptomatic peripheral end-neuromas

NEUROCAP[®], synthetic and bioresorbable, is the only approved nerve capping device targeted for surgical management of symptomatic end-neuromas. NEUROCAP[®] is intended to protect a peripheral nerve-end and to reduce the development of a symptomatic neuroma.

Interim animal data support the mechanism of action of NEUROCAP[®] as an effective nerve-end barrier

In our ongoing animal study, rat sciatic nerve defects were treated by capping the nerve-end with NEUROCAP[®] compared to a control group that was treated with the cut and bury technique. The effect of NEUROCAP[®] on the nerve-end and its surrounding tissues will be histologically assessed at 3, 6 and 12 months after implantation on myelin content and axon evaluation, neuroma formation, nerve outgrowth, abnormal healing of the nerve, and any abnormal response in the surrounding nerve tissues.



Preliminary 6-month histological results show NEUROCAP[®] is safe to use and effective in providing the expected barrier function for the nerve-end, while sprouting and tethering is inhibited.

Schematic representation of the obtained results



NEUROCAP[®]



Histological evaluation at 3 and 6 months after implantation (10x total magnification)



NEUROCAP[®] effectively provides the expected barrier function for the nerve-end, while sprouting is inhibited. The right picture is a graphical illustration of the histological assessment. **Dark Blue**: outline of NEUROCAP[®] implant at 3 and 6 months post-op (rat-sciatic nerve); note the extended device degradation at 6 months. **Turquoise Blue**: collaborated inner part of the NEUROCAP[®] (no nerve tissue present), marginally observed at 6 months due to further device degradation **Orange**: Fixating sutures; **Magenta Red**: outline of nerve-end

CONTROL GROUP

Histological evaluation at 3 and 6 months after implantation (10x total magnification)



The control group, where the barrier is missing, allows for nerve-end outgrowth. The right picture is a graphical illustration of the histological assessment at 3 and 6 months post-op. **Dark Blue**: outline of nerve-end; **Magenta Red**: outline of newly formed neuroma and fibrotic tissue; **Orange**: Fixating suture

Polyganics - Rozenburglaan 15A, 9727 DL Groningen, The Netherlands
T +31 (0)50 588 65 88 - info@polyganics.com - www.polyganics.com

The information presented in this brochure is intended to inform on the product. Always refer to the package insert, product label and/or user instructions before using this product. NEUROCAP is a registered trademark of and manufactured by Polyganics B.V., The Netherlands. This biocompatible device is composed of proprietary compositions of copolyester poly-lactide-caprolactone, PLCL.

POLYGANICS

TRANSFORMING PATIENT RECOVERY