Instructions for use - GB

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications

STERILE

Sterilized with ethylene oxide gas. For single use only. Do not autoclave.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

STORAGE

- Store in a dark, dry place between -18°C and 8°C (39°F).
- Use the device prior to the 'use by' specified on the package.

Description

The Vivosorb® sheet is composed of a bioresorbable copolyester poly(DL-lactide-ε-caprolactone) sheet. The Vivosorb® sheet can be applied to temporarily separate opposing soft tissues.

The Vivosorb® sheet elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the material by fibrous tissue. Degradation of the sheet occurs through hydrolysis leading to gradual reduction of molecular weight. The sheet remains its initial mechanical properties up to at least 8-10 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and ω-hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. The material is resorbed within approximately 16 months.

The Vivosorb® size and thickness are indicated on the label. Vivosorb® is packed in a Tyvek pouch protected by a Tyvek sheet. The Tyvek pouch is sealed in an aluminium pouch including silica gel, and then packed in a carton box with an IFU.

Vivosorb® is indicated for single-use.

Indications

Vivosorb® is indicated for use as a temporary protective sheet to separate opposing soft tissues, thereby minimizing tissue attachment.

Contraindications

There are no known contraindications.

Warnings

- The Vivosorb® sheet is for single use only. Do not resterilize or re-use. Structural integrity and/or function may be impaired through cleaning, resterilization, or re-use and may cause adverse patient reactions. Accordingly, Polyganics will not be responsible for any direct or consequential damages or expenses resulting from re-use of (or any part of) the Vivosorb® sheet;
- Vivosorb should not be used in cardiac surgery or used in direct contact with the heart, the central circulatory system or the central nervous system;
- Sterile unless package has been opened or damaged. Discard open unused devices;
- The Vivosorb® sheet should only be used by physicians who are trained in surgery. Accordingly, Polyganics will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel. The physician should consult recent literature on current medical practice on surgery.
- The Vivosorb® sheet is not intended to bridge tissue gaps.
- Tissue support may be suboptimal in elderly, malnourished or debilitated patients or in patients suffering from cancer, anaemia, obesity, diabetes,

infection or other conditions which may delay wound healing, infected wounds, or mild tissue inflammatory response characteristic of foreign body response.

Precautions

- Do not expose the sheet to organic solvents (e.g. chloroform, acetone);
- Always fixate the sheet with non absorbable sutures or absorbable sutures keeping their strength up to at least 12 weeks
- Wrinkling or folding may cause sticking of the sheet to itself.
- Ensure sufficient healthy soft tissue is available to cover VIVOSORB in order to avoid protrusion or wound dehiscence.

Adverse Events

Procedures requiring implantation of a Vivosorb® sheet should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during the procedure.

Adverse events associated with the use of a Vivosorb® sheet may include but are not limited to:

- Failure to provide adequate support for appropriate healing;
- Transitory local irritation;
- Infection;
- Allergy;
- Delayed wound healing.

Opening of the package


Open the aluminium pouch. **Attention:** the outside of white/transparent pouch is not sterile.

The pouch is opened in such a way that the device remains sterile. By claspings the sheet gently by forceps, it can be taken from the tyvek protective sheet.

Suturing Technique

1. Fixate the Vivosorb® sheet to the surrounding tissue (using standard suturing techniques). Ensure that there is sufficient soft tissue to cover the device.
2. Or fixate simultaneously with closing the wound. The later option makes that the sheet is an integral part of the closed tissue and sutures.
3. The sutures must be at least 2mm from the edges of the sheet.
4. It must be ensured that the sheet does not act as a bridge to close the gap.
5. Dispose contaminated implantation and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

USED SYMBOLS

-  = Consult instructions for use.
-  = Caution
-  = Do not use if package is damaged
-  = Manufacturer
-  = Use by
-  = Catalog number
-  = Do not reuse
-  = Do not re-sterilize
-  = Lot number
-  = Sterilized with ethylene oxide gas
-  = Keep dry
-  = Keep away from sunlight
-  = Storage temperature limit

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Polyganics product(s) described in this publication. Under no circumstances shall Polyganics be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Polyganics to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Polyganics printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Polyganics will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Manufacturer:
 Polyganics Tel: +31 50 588 6588
 Blauwborgje 32 www.polyganics.com
 9747 AC Groningen
 The Netherlands

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